

Aspire VIDA Folding Wheelchair – Transit

SKU: MWS449850 (Orange)

SKU: MWS449855 (Green)

\$675



Features and Benefits:

- The Aspire Vida is a super lightweight transport wheelchair ideal for quick trips to the shops, appointments and catching up with friends.
- Ergonomic backrest and seat shape follow the contours of the user to maximise comfort and postural support.
- The 135kg safe working load caters to a broad range of users.
- Lightweight frame design ensures easy lifting in and out of cars, ideal for carers.
- Contoured and padded seating improves comfort and postural support.
- Attendant operated hand brakes with parking lock.
- Simple and intuitive folding for easy transportation.
- Removable leg rests and folding back for compact storage.
- Rear storage pocket for essentials.
- Available as attendant propelled or self-propelled.
- Anterior Seat Height - 495 mm
- Posterior Seat Height - 445 mm

Specifications:

Total Length	995mm	Colour	Orange, Green
Total Width	655mm	Seat Depth	420mm
Total Height	930mm	Seat Width	450mm
SWL	135kg	Frame Weight	8.6kg
Unit Weight	9.6kg	Backrest Height	430mm