

Shower Stool With Padded Seat

SKU: KA222ZAP

\$130



Features and Benefits:

- Designed for users preferring additional seating support while showering .
- White powder coated 22mm steel tube .
- Adjustable height and angle to suit individual users and uneven floors .
- Padded seat clips on for easy assembly, change and cleaning .
- Anti-slip rubber tips on each leg for stability .
- Independently tested to AS/NZS 3973:2009 .
- Suitable for institutional and domestic use .
- Available in 45cm and 50cm seat widths .

Specifications:

Code	KA222ZAP	Max User Weight	110kg
Overall Width	510mm	Overall Depth	380mm
Seat Width	340mm	Seat Depth	290mm
Min Seat Height	510mm	Max. Seat Height	630mm